

DNA extraction from sterivex filter
(Staley Lab 2008)

Modified from:

Vetriani et al., (2003) Applied and Environmental Microbiology 69: 6481-6488.

Solutions needed:

Solution 1

50 mM glucose

10 mM EDTA

25 mM Tris-Cl (pH 8)

Lysozyme Solution (Remake daily)

0.8 mg/mL lysozyme in Solution 1

Proteinase K Solution (Remake daily)

10 mg/mL in TlowE

- 1.) Heat close filter bottom. Parafilm bottom.
- 2.) Pipet 1 mL of Solution 1 onto filter.
- 3.) Freeze thaw sample between dry ice/Ethanol slurry and 55°C water bath 8-10 times. Transfer quickly.
- 4.) Add 400 μ L Solution 1, 200 μ L Lysozyme solution, 100 μ L 500 mM EDTA, 20 μ L proteinase K. Parafilm top.
- 5.) Turn for 30 minutes at 30°C.
- 6.) Take solution out of filter case with syringe. Split between 2 or 3 microcentrifuge tubes.
- 7.) Add 150 μ L of 10% SDS followed quickly by phenol chloroform (pH 7) (1:1 volume with sample). Invert 10 times.
- 8.) Spin in centrifuge at max speed for 3 minutes.
- 9.) Remove top phase with pipet. Add to new tube.
- 10.) Add 1:1 volume of phenol chloroform. Invert 10 times. Spin in centrifuge at max speed for 3 minutes. Remove top phase with pipet. Add to 2 mL tube.
- 11.) Precipitate DNA: add 50 μ L of 3.0 M sodium acetate and 2X volume of 100% ethanol.
- 12.) Put tubes in -15°C freezer for at least 3 hours. (Overnight is ok.)
- 13.) Spin in centrifuge at max speed for 15 minutes.
- 14.) Wash pellet with 100 μ L 70% ethanol. Centrifuge 5 minutes on max.
- 15.) Dry in speed-vac.
- 16.) Resuspend in 100 μ L TlowE for several hours in the fridge.
- 17.) Clean DNA with Qiagen Qiaquick PCR purification kit.
- 18.) Resuspend in 50 μ L TlowE.