Salt DNA extraction (C. Lovejoy)

Note: The numbering is our protocol for 0.22 μm Sterivex (Millipore) filters, The step b is for cases when we used 47 mm polycarbonate filters, .

Step:

0. Put 70% ethanol in the freezer. You need \pm 6 ml per sample.
1. Thaw your sample (RT or 37°C).
 - 1b. Transfer filters and lysis buffer in a 15 ml tube.
2. Add 50 μl of lysozyme, incubate 45 min at 37°C, with rotation.
3. Add 50 μl of proteinase K and 200 μl of SDS 10%, incubate 1h at 55°C, with rotation.
4. Transfer sterivex content in a new 15ml tube and add 1ml of lysis buffer in the sterivex. Incubate for 15 min at 55°C with rotation and then transfer to the 15 ml tube.
 - 4b. For 47 mm filters, add 1 ml of lysis buffer in the 15 ml tubes.
5. Add 1.9 ml of NaCl 6M. Vortex 1 min. Centrifuge at maximum speed for 20 minutes.
6. Transfer flow through in a new 15 ml tube (At this step you there you can often see a white fuzzy form, you can centrifuge longer. Be careful to not transfer too much of that in the new tube.).
7. Add 5 ml of cold 70% ethanol, mix by inverting the tubes and put at -20 (or on the ice) for at least 10 minutes.
8. Transfer 2 ml of this solution in 2 ml tubes (you can use two 2 ml tubes per sample to go faster).
9. Spin 10 minutes at maximum speed and discard flow through.
10. Using the same 2 ml tube, repeat steps 8-9 until all the solution has been centrifuge.
11. Quick spin to remove remaining ethanol. Be careful to not remove the pellet, it can be invisible.
12. Add 200 μl of cold ethanol. Spin 5 minutes at maximum speed. Remove ethanol. Let pellets dry (no alcohol trace, but do not let dry for too long).
13. Add 100ml of TE, stand for 1-5 min. Dissolve by up and down and vortex. Transfer in a clean tube.

14. Make aliquot, check on gel and quantify. Keep DNA samples at -20 or -80.

Solutions needed:

70% ethanol, for 500 ml
mix 368,4 ml of 95% ethanol with 131,6 ml of milliQ water.

Lysozyme (SIGMA Cat# L-7651) 2 mg lysozyme in 50 µl lysis buffer. Prepare fresh.

Proteinase k (SIGMA Cat# P-2308) 0.4 mg proteinase K in 50 µl lysis buffer. Prepare fresh.

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