

Seawater nucleic acid extraction protocol vs. 1.4

1. Add 1.6 ml of SET lysis buffer directly into top of Sterivex using a 2.5-ml syringe with a 25G 5/8" needle.
2. Add 180 μ l of fresh lysozyme.
3. Seal large end with a cut 1-ml syringe tip (containing plunger fragment).
4. Cover narrow end with blue tac.
5. Incubate at 37°C for 30 min with rotation in a hybrid oven.
6. Add 200 μ l of SDS.
7. Add 55 μ l 20 mg/ml fresh proteinase K.
8. Incubate at 55°C for 2 h with rotation in a hybrid oven.
9. Withdraw lysate into a 5 ml syringe.
10. Add 1 ml SET buffer to Sterivex and rotate to rinse.
11. Withdraw rinse buffer into same 5 ml syringe.
12. Add lysate to 15-ml Phase-Lock tube containing 2 ml of phenol:chloroform:isoamyl alcohol (25:24:1), pH 8 (in Maxtract gel lock system, QIAGEN). Shake gently until mixed.
13. Centrifuge @ 1,500 x g for 5 min.
14. Add an additional 2 ml of phenol:chloroform:isoamyl alcohol (25:24:1). Shake gently until mixed.
15. Centrifuge @ 1,500 x g for 5 min.
16. Add 2 ml of chloroform:isoamyl alcohol (24:1). Shake gently until mixed.
17. Centrifuge @ 1,500 x g for 5 min.
18. Decant aqueous phase to a sterile and DEPC-treated (if RNA needed) JA-20 centrifuge tube.
- 18b. Add 5 μ l of glycogen (Roche).
19. Add 0.5V of 7.5 M NH₄acetate. Mix briefly.
20. Add 2.5V of pure ethanol.
21. Mix and leave @ -20°C for > 1 h (overnight is fine).

22. Centrifuge in a JA-20 rotor @ 20,000 rpm for 30 minutes @ 4°C.
23. Decant ethanol.
24. Add 2 ml 80% ethanol and rinse tube.
25. Centrifuge in a JA-20 rotor @ 20,000 rpm for 20 minutes @ 4°C.
26. Decant ethanol.
27. Add 2 ml 80% ethanol and rinse tube.
28. Centrifuge in a JA-20 rotor @ 20,000 rpm for 20 minutes @ 4°C.
29. Decant ethanol and leave inverted.
30. Dry inverted tube for 15 minutes in fume hood (provides air flow).
31. Suspend invisible pellet in 200 µl DEPC-treated sterile water. Leave on ice for approximately 1 h with frequent finger-tapping to rinse tube walls.
33. Store @ -20°C (DNA) or -80°C (RNA).

SET buffer

40 mM EDTA

50 mM Tris-HCl, pH 9

0.75 M sucrose

autoclave and store @ RT

Lysozyme

990 µl sterile water

9 mg lysozyme

9 µl 1M Tris-HCl, pH 8

store on ice

SDS

10% SDS

Proteinase K

950 μ l sterile water

50 μ l 1M Tris-HCl, pH 8

20 mg proteinaseK

7.5 M NH₄acetate

According to molecular cloning manual