

DNA extraction from Cultured Algae (Robert Anderson and HwanSu Yoon, Bigelow Laboratory)

We prepared two tubes of sample (SAB_0001_2009_01_28 and SAB_0002_2009_01_28).

Samples were prepared as follow;

- 1) We removed 100 uL from each culture into 1.5 mL eppendorf tubes (total 18 tubes from ca 300 cultures).
- 2) These tubes were spun down at high speed to pellet the bacteria.
- 3) We extracted the DNA from each 18 tubes separately using Qiagen DNeasy Plant mini kit.
- 4) We did PCR using bacterial 16S rDNA specific primers. Most of reactions have positive band (see below gel picture).
- 5) Then combined extracted DNA from 18 samples.
- 6) One half of DNA was purified again using Qiagen QIAamp DNA mini kit (SAB_0001_2009_01_28). DNA concentration was 26.4 ng/ul. Total 90 ul was left on the tube (2,376 ug, total).
- 7) Concentration of rest of DNA (500 ul) was 5 ug/ul (SAB_0001_2009_01_28).
- 8) Two DNAs were dried after a normal EtOH precipitation method (with 3M NaOAc).